

# RNA extraction for BLiMMP project

*Written by Benjamin D. Peterson*

## Solutions/supplies needed

- Lysis matrix A in MP Bio bead-beating tubes
- TRIzol
- Internal standard, at ~0.75-1.25 ng/ $\mu$ l
- Nuclease-free 1.5 ml Eppendorf tubes (4 per extraction)
- Chloroform
- 100% isopropanol
- 100% ethanol, cold
- 70% ethanol, cold
- 3M sodium acetate, pH = 5.5
- Turbo DNase kit

### **Prep work before starting extractions**

- 1 Cool centrifuge to 4°C.
- 2 Set up and label four 1.5 ml Eppendorf tubes for each extraction
- 3 Clean everything with RNase Zap.

### **Physical lysis and Trizol extraction**

- 1 Sterilize and RNase Zap the PVC cutter, a razor blade, and two tweezers.
- 2 Cut open Sterivex filter cartridge with the PVC cutter and remove the filter using the razor and the tweezers.
- 3 Add filter to beadbeating tube with lysis matrix A in it. Immediately add 1 ml of Trizol and bead-beat for 1 minute.
  - a. While the sample is bead-beating, wipe down and RNase Zap the tweezers, razor, and PVC cutter
- 4 Store sample in fridge.
- 5 Repeat for all the samples that need processing on this day.
- 6 Now process all the samples together. Add 12 ng of internal standard to each sample. For the internal standard prepped on 2022-03-29, that's 14 µL of the 500X diluted internal standard. Invert repeatedly for 1 min.
- 7 Centrifuge for 5 minutes at 16,100 x g at 4°C .
- 8 Transfer 800 µl of supernatant to new 1.5 ml Eppendorf tube.
- 9 Add 300 µL of the chloroform to Trizol/sample mix. Invert repeatedly for 90 seconds, then incubate at room temperature for 3 minutes.
- 10 Centrifuge for 15 minutes at 13,200 x g at 4°C.
- 11 Transfer 500 µL of the aqueous phase to a new 1.5 ml Eppendorf minicentrifuge tube.

### **Purification by isopropanol and ethanol precipitation**

- 12 Add 500 µL of 100% isopropanol to sample. Invert several times to mix.
- 13 Incubate at room temperature for 20 minutes.
- 14 Spin down at 16,100 x g for 25 minutes at 4°C.
- 15 Carefully remove the supernatant by pipetting it off.
- 16 Resuspend pellet in 100 µl of DEPC-treated H<sub>2</sub>O.
- 17 Add 10 µl of 3 M NaOAc at pH 5.5. Mix well.
- 18 Add 250 µl of cold 100% EtOH. Mix well, store at -20°C for 20 minutes.
- 19 Centrifuge for 20 min at 16100xg at 4°C.
- 20 Add 200 µl of 70% EtOH, and vortex to mix.
- 21 Spin down at 16,100 x g for 20 minutes at 4°C.
- 22 Carefully pipette off the supernatant and allow RNA to dry on the countertop.
- 23 Resuspend in 50 µl of nuclease-free water.

### **DNA removal with Turbo DNase**

- 24 Add 5 µl DNase Buffer from Turbo DNA kit. Mix gently
- 25 Add 1.5 µl of TURBO DNase. Mix gently. Incubate at 37 °C for 25 minutes.
- 26 Add 5µl of DNase Inactivation Reagent (mix reagent well first). Incubate 5 min at room temp, mixing periodically.
- 27 Centrifuge for 90 seconds. Transfer supernatant to new tube, take ~52 µl. Be sure not to get any of the inactivation reagent.

### **Purification by ethanol precipitation**

- 28 Add 148 µl of nuclease free H<sub>2</sub>O to sample, mix well
- 29 Add 20 µl of 3 M sodium acetate (pH 5.5) to the sample and mix well.
- 30 Add 500 µl of cold 100% ethanol to sample, mix well.
- 31 Store at -20°C for 20 minutes.
- 32 Centrifuge sample for 20 minutes at 16,100 RCF and at 4°C.
- 33 Pipet off supernatant.
- 34 Add 350 µl of cold 70% ethanol. Vortex briefly.
- 35 Centrifuge sample for 20 minutes at 16,100 RCF and at 4°C.
- 36 Pipet off supernatant, allow RNA to dry on bench top.
- 37 Resuspend in 54 µl of nuclease-free water.
- 38 Remove 4 µl aliquot for Qubit (both DNA and RNA) and Nanodrop, immediately store the remainder at -80°C.