

Prepare the internal standard for RNA extractions.

Solutions/supplies needed

- pFN18a plasmid
- BamHI
- Buffer E (for BamHI)
- BSA (10ug/uL)
- Water bath at 37°C
- QIAquick PCR purification
- Mung bean nuclease
- Mung bean nuclease buffer
- Glycerol (100%)
- Riboprobe System T7
 - Transcription Optimized 5X buffer
 - 100mM DTT
 - Recombinant RNasin Ribonuclease Inhibitor (100u total)
 - Nucleotide mix
 - RNA polymerase
 - DEPC water
- Qiagen RNeasy Mini kit
 - Buffer RLT
 - 100% ethanol
 - RNeasy spin column
 - Buffer RW1
 - DNase i
 - RPE buffer

Protocol

1. Linearize the pFN18a plasmid with PmeI

- Mix
 - 7 uL H₂O
 - 2 uL CutSmart Buffer (10x)
 - 10 uL pFN18A plasmid (100 ng/uL)
 - 1 uL PmeI
- Incubate at 37°C for 90 minutes, then inactivate at 65°C for 20 minutes

2. Extract and purify DNA by phenol:chloroform extraction and ethanol precipitation

- Move 20μl of digestion mix to a 1.5ml Eppendorf microcentrifuge tube
- Add 50μl of phenol:chloroform:IAA (pH~8)
- Invert or shake sample for 2 minutes
- Centrifuge at max speed for 5 minutes
- Remove upper aqueous layer into a new 1.5ml microcentrifuge tube

- Add 2 μ l of filter-sterilized 3M NaAc (pH 5.5) and mix
- Add 200 μ l of 100% EtOH and invert to mix
- Precipitate DNA for 2hr at -20°C.
- Spin down the sample for 30 minutes at top speed
- Pipette off the supernatant
- Rinse with 300 μ l of 70% EtOH, invert to mix
- Spin down for 25 minutes at max speed
- Pipette off the supernatant
- Allow DNA to dry
- Resuspend in 20 μ l of NF-H₂O

3. Run in vitro transcription kit (Riboprobe System T7)

- Recipe:
 - 20 μ l Transcription Optimized 5X buffer
 - 10 μ l 100mM DTT
 - 2.5 μ l Recombinant RNasin Ribonuclease Inhibitor (100u total)
 - 20 μ l Nucleotide mix
 - To prep this, mix equal parts of the four different 10mM rNTPs stocks
 - 2 μ l Linearized template DNA (1-2.5mg/ml)
 - 2 μ l T7 RNA polymerase (40u total)
 - 43.5 DEPC-treated water (bringing up to 100 μ l)
- Incubate mix for 1-2 hours at 37-40°C

4. DNA removal and RNA purification (Qiagen RNeasy Mini kit)

- Prep solutions:
 - Add 4 volumes of 100% EtOH to Buffer RPE (only needs to be done once)
- Add 350 μ l Buffer RLT and mix
- Add 250 μ l 100% ethanol and mix well
- Transfer sample to RNeasy spin column
- Centrifuge 15s at >10,000rpm then discard flow-through
- Add 350 μ l of Buffer RW1 to column, spin at >10,000rpm for 15s
- Add 10 μ l DNase I to 70 μ l Buffer RDD, mix and centrifuge
- Add 80 μ l DNase/RDD mix to RNeasy spin column.
- Incubate at room temperature for 15 min.
- Add 350 μ l RW1 buffer to RNeasy spin column, spin at >10,000rpm for 15s, discard flow-through
- Add 500 μ l Buffer RPE to RNeasy spin column, spin at >10,000rpm for 15s, discard flow-through
- Spin RNeasy spin column at >10,000rpm for 1min
- Place RNeasy spin column in new 1.5ml Eppendorf centrifuge tube
- Add 50 μ l of DEPC-treated water directly to RNeasy spin column, spin at >10,000rpm for 1min

9. Qubit quantification + run in gel